

OUR NEW EQUIPMENT



Tableting press simulator **STYL'ONE EVOLUTION** Tablet press for R&D and clinical batches manufacturing

Formulation development

- . Impact of formulation on dissolution
- . Impact of tooling characteristics

Screening/Study of raw materials and blends

- . Compression cycle studies
- . Pressure/Hardness relationship
- . Compaction profiles

Process development & trouble shooting investigation

- . Determination of compression parameters & critical parameters
- . Simulation of scale-up parameters & industrial tableting machines

Manufacturing of technical & small clinical batches

- . From 150 g to 5 kg



Evoluteive on demand

(multi-layer, roller compaction simulator, Tab-in-Tab)



For further information, please get in touch with Olivier Cadiou.

olivier.cadiou@synerlab.com

Our core business:

- Pharmaceutical and Analytical Development of innovative drugs, including manufacturing and control of clinical batches.

Our certifications:

- Certificate of compliance with current Good Manufacturing Practices (cGMP) for medicinal products & investigational medicinal products for human & veterinary use.
- Accredited to research tax credit and innovation tax credit.